

<b>AUTOCLAVE MAINTENANCE LOG</b>			
<b>WEEK</b>	<b>DATE</b>	<b>BIOLOGICAL INDICATOR RESULT (POSITIVE/NEGATIVE?)</b>	<b>INITIALS</b>
<b>WEEK 1</b>			
<b>WEEK 2</b>			
<b>WEEK 3</b>			
<b>WEEK 4</b>			
<b>WEEK 5</b>			

<b>WEEKLY</b>		
<b>DATE</b>	<b>ACTION</b>	<b>INITIALS</b>
	CLEAN CHAMBER AND TRAYS, DRAIN WATER AND REFILL THE RESIVOR WITH DISTILLED OR DEMINERALIZED WATER.	

<b>MONTHLY</b>		
<b>DATE</b>	<b>ACTION</b>	<b>INITIALS</b>
	FLUSH THE SYSTEM WITH SPEED CLEAN STERILIZER CLEANER. (SEE INSTALLATION & OPERATION MANUAL)	

ANY ISSUES WITH BIOLOGICAL INDICATOR HAVING INADEQUATE RESULTS OR ANY OTHER PROBLEMS WITH THE AUTOCLAVE PLEASE SEE CLINICAL COORDINATOR

MAINTENANCE MUST BE DONE IN ORDER TO STAY IN COMPLIANCE